



SYNTERVENTION™

SWIPER™

medical device foam wiper

Qty: 1

Size: 3" x 3" x 0.25"

Instructions for Use:

1. Using aseptic technique, peel open package and transfer Swiper to sterile field/personnel.
2. Hydrate Swiper by compressing, submerging and releasing in sterile saline/water.
3. With light pressure, wipe medical device with hydrated Swiper.

For more information or to place an order visit

www.syntervention.com

REF 5001

LOT 123456

2010-05

STERILE | R



Made in Thailand for
Syntervention, Inc.
1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804
888-505-8802 main
888-853-6559 fax

00024806 Rev. B



(01) 0 0858303 00300 2 (17) 100531 (10) 123456



SYNTERVENTION™

SWIPER™

medical device foam wiper

Qty: 2

Size: 3" x 3" x 0.25"

Instructions for Use:

1. Using aseptic technique, peel open package and transfer Swiper to sterile field/personnel.
2. Hydrate Swiper by compressing, submerging and releasing in sterile saline/water.
3. With light pressure, wipe medical device with hydrated Swiper.

For more information or to place an order visit

www.syntervention.com

REF 5002

LOT 123456

2010-05

STERILE R



Made in Thailand for
Syntervention, Inc.
1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804
888-505-8802 main
888-853-6559 fax

00024806 Rev. B



(01) 0 0858303 00301 9 (17) 100531 (10) 123456



SYNTERVENTION™

SWIPER™

medical device foam wiper

Qty: 4

Size: 3" x 3" x 0.25"

Instructions for Use:

1. Using aseptic technique, peel open package and transfer Swiper to sterile field/personnel.
2. Hydrate Swiper by compressing, submerging and releasing in sterile saline/water.
3. With light pressure, wipe medical device with hydrated Swiper.

For more information or to place an order visit

www.syntervention.com

REF 5004

LOT 123456

2010-05

STERILE R



Made in Thailand for
Syntervention, Inc.
1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804
888-505-8802 main
888-853-6559 fax

00024806 Rev. B



(01) 0 0858303 00302 6 (17) 100531 (10) 123456



SYNTERVENTION™

SWIPER™

medical device foam wiper

Qty: 100, 1/PKG

Size: 3" x 3" x 0.25"

Made in Thailand for
Syntervention, Inc.
1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804
888-505-8802 main
888-853-6559 fax

LB3007 Rev.0A

REF 5001

LOT 123456

2010-05

STERILE R



For more information or to place an order visit
www.syntervention.com



(01) 1 0858303 00300 9 (17) 100531 (10) 123456



SYNTERVENTION™

SWIPER™

medical device foam wiper

Qty: 100, 2/PKG

Size: 3" x 3" x 0.25"

Made in Thailand for
Syntervention, Inc.
1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804
888-505-8802 main
888-853-6559 fax

LB3007 Rev.0A

REF 5002

LOT 123456

2010-05

STERILE R



For more information or to place an order visit
www.syntervention.com



(01) 1 0858303 00301 6 (17) 100531 (10) 123456



SYNTERVENTION™

SWIPER™

medical device foam wiper

Qty: 100, 4/PKG

Size: 3" x 3" x 0.25"

Made in Thailand for
Syntervention, Inc.
1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804
888-505-8802 main
888-853-6559 fax

REF 5004

LOT 123456

2010-05

STERILE R



L23007 Rev.0A

For more information or to place an order visit
www.syntervention.com



(01) 1 0858303 00302 3 (17) 100531 (10) 123456