

Swasher Ultra-Low-Lint Surgical Towels

Syntervention
(888) 853-6559
www.syntervention.com

KEY FEATURES

- 99% less lint than standard woven towels
- Excellent fluid capacity for uncompromised performance
- Lightweight to reduce medical waste
- Latex- and DEHP-free

Swasher is an ultra-low-lint and low-particulate surgical towel from Syntervention (Rocky Mount, NC) designed to reduce lint and particulate contamination and minimize the transfer of foreign materials in inter-

ventional and surgical procedures. Swasher weighs less than a regular towel, with the potential to reduce the cost of medical waste disposal. It is biocompatible, latex- and DEHP-free, and available for both sterile and nonsterile applications. The product is soft and drapable, yet strong, with a "grab tab" fold for easy opening.

