

GUIDE-KLIP™




MEDICAL DEVICE RETAINING TOOL

APPLICATIONS

- Diagnostic and interventional procedures
- Operating Room
- ICU/CCU
- Emergency Room

FEATURES

- Adhesive Pull Tab
- Water resistant
- Design allows for devices to move freely thru fenestrations



The GUIDE-KLIP is a sterile, water resistant retention tool designed to reduce accidental contamination by holding and organizing guidewires, catheters, and cables in the patient care setting.



SYNTAX®
Bringing Synergy to Intervention

GUIDE-KLIP Ordering Information is also available on our website.



"Not For Human Use"

Part Number	Description	Carton Quantity
5100	GUIDE-KLIP, 3 per pack	100 units/carton

Remit To Syntervention, Inc., 1444 Jeffreys Road, Suite 133, Rocky Mount, NC 27804

Contact Us Phone: 888.505.8802
 Fax: 888.853.6559
 Email: customerservice@syntervention.com

Tax ID 45-3609564

Shipment Items are shipped ground unless otherwise specified.

Product Storage Store in a cool, dry place. Avoid exposure to moisture and excessive heat.

Product Sterilization The GUIDE-KLIP is sterilized with EO.

Shelf Life See package label for expiration date.

Latex and DEHP Status GUIDE-KLIP is not made with natural rubber latex and does not contain DEHP.

Returned Goods Policy A Returned Goods Authorization (RGA) number is required for all returns and replacements. Authorization can be obtained by calling Syntervention Customer Service at 888.505.8802.



Syntervention, Inc.
 1444 Jeffreys Road, Suite 133
 Rocky Mount, NC 27804
 Customer Service: 888.505.8802
 Fax: 888.853.6559
www.syntervention.com