


SYNTERVENTION

SWIPER[®]

Lint-Free Medical Device Foam Wiper



The SWIPER Lint-Free Medical Device Foam Wiper

The SWIPER is a sterile, biocompatible, cellulose-free foam wiper designed for the removal of blood, contrast and other contaminants from sterile instruments and medical devices.

Improved Efficiency

- Absorbs 20x its weight in saline solution.
- Transfers fluid evenly to medical device surface.
- Durable and non-shredding.

Less Waste

- Two (2) SWIPERS are the typical amount of product used during a procedure.
- Use our calculators to estimate how much your facility could save. Visit syntervention.com/calculators.

Effortless Handling


- Designed for easy, repetitive manipulation.
- Conformable yet maintains shape when wet.
- Fits comfortably into any hand size.

SWIPER has 93% less lint than gauze and TELFA.*

Product	≥ 10 μm	≥ 25 μm	≥ 50 μm	≥ 100 μm
Cotton Gauze	29,809	937	8	0
TELFA	27,884	2,401	111	15
SWIPER	2,090	60	2	0

*Independent LPC testing. Data on file at Syntervention, Inc.

Lose the Lint. Switch to the SWIPER.





SWIPER Order Information

Number	Description	Size	Case Quantity
5001	SWIPER, 1 per pack	3 x 3 x 0.25	100 units/case
5002	SWIPER, 2 per pack	3 x 3 x 0.25	100 units/case

Remit To	Syntervention, Inc., 1444 Jeffreys Road, Suite 133 Rocky Mount, NC 27804
Contact Us	Phone: 888.505.8802 Fax: 888.853.6559 Email: info@syntervention.com
Tax ID	45-3609564
Payment Terms	Net 30 days.
Shipment	Items are shipped ground unless otherwise specified.
Product Storage	Store in a cool, dry place. Avoid exposure to moisture and excessive heat.
Product Sterilization	Sterilized with EO or gamma radiation. Non-sterile products available upon request.
Shelf Life	See package label for expiration date. 1-3 years.
Latex and DEHP Status	Latex, BPA and DEHP free.
Returned Goods Policy	A Returned Goods Authorization (RGA) number is required for all returns and replacements. Authorization can be obtained by calling Syntervention Customer Service at 888.505.8802.

To view calculators visit www.syntervention.com/calculators.

To view the Learning Center Video Library visit www.syntervention.com/video-library.

Syntervention and SWIPER, are trademarks or registered trademarks of Syntervention, Inc. TELFA is a trademark of Covidien AG or an affiliate.



SYNTERVENTION

Syntervention, Inc.
1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804
Customer Service: 888.505.8802
Fax: 888.853.6559
www.syntervention.com