


Sterile Product – OPERATING ROOM PROCEDURES

SWICKER Part Number	Dimensions (in)	Thickness (in)	Firmness	Radiopaque dots per sponge	Quantity	Packs Per Case
4008-203	1.5 x 3	.125	Soft	5	2/pack	300
4008-303	1.5 x 3	.25	Medium	5	2/pack	300
4008-403	1.5 x 3	.50	Firm	5	2/pack	200
4008-205	3 x 3	.125	Soft	6	2/pack	300
4008-305	3 x 3	.25	Medium	6	2/pack	200
4008-405	3 x 3	.50	Firm	6	2/pack	100
4008-207	5 x 5	.125	Soft	10	2/pack	150
4008-307	5 x 5	.25	Medium	10	2/pack	100
4008-407	5 x 5	.50	Firm	10	2/pack	50
4008-209	9 x 9	.125	Soft	18	2/pack	75
4008-309	9 x 9	.25	Medium	18	2/pack	50
4007-409	9 x 9	.50	Firm	18	1/pack	50
4011	DIGIT	.25	Medium	3	1/pack	100
4013	MITT	.25	Medium	4	1/pack	100



Sterile Product – LAPAROSCOPIC AND ROBOTIC PROCEDURES

SWICKER Part Number	Dimensions (in)	Thickness (in)	Firmness	Radiopaque dots per sponge	Trocar size (mm)	Quantity	Packs Per Case
4008-203	1.5 x 3	.125	Soft	5	≥ 5	2/pack	300
4008-303	1.5 x 3	.25	Medium	5	≥ 10	2/pack	300
4008-403	1.5 x 3	.50	Firm	5	≥ 12	2/pack	200

SWICKER®

Radiopaque Surgical Foam

Remit To: Syntervention, Inc., 1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804

Contact Us: Phone: 888.505.8802
Fax: 888.853.6559
Email: info@syntervention.com

Tax ID: 45-3609564

Payment Terms: Net 30 days.

Shipment: Items are shipped ground unless otherwise specified.

Product Storage: Store in a cool, dry place. Avoid exposure to moisture and excessive heat.

Product Sterilization: The SWICKER is sterilized with EO.

Product Shelf Life: See package for Expiration Date.

Hazardous Material Status: Latex, BPA and DEHP free.

California Requirements: Meets Proposition 65 standards.

SWICKER Non- Sterile: The SWICKER is available in non-sterile version.

Returned Goods Policy: A Returned Goods Authorization (RGA) number is required for all returns and replacements. Authorization can be obtained by calling Syntervention Customer Service at 888.505.8802.



SPONGE IFU



DIGIT IFU



MITT IFU



SYNTERVENTION

Syntervention, Inc.
1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804
Customer Service: 888.505.8802
Fax: 888.853.6559
www.syntervention.com

To view calculators visit www.syntervention.com/calculators.

To view the Learning Center Video Library visit
www.syntervention.com/video-library.

Syntervention and SWICKER, are trademarks or registered trademarks of Syntervention, Inc.